

# Implant Prosthetic Components – Instructions for Use

IMPORTANT INFORMATION — PLEASE READ

## Caution

U.S. federal law restricts this device to sale by, or on the order of, a licensed dentist or physician.

## General Information

PREAT prosthetic components are for use by qualified, licensed clinicians and laboratory technicians fully trained in their application.

## Disclaimer of Liability

The guidelines presented herein are not adequate to allow inexperienced clinicians to administer prosthetic dentistry, and are not intended to substitute for formal clinical or laboratory training.

These devices should only be used by individuals with training and experience specific to their clinically accepted application.

PREAT is not liable for damages resulting from treatment outside of our control. The responsibility rests with the provider.

## Prosthetic Considerations

Following successful placement of dental implants, verify primary stability and appropriate occlusal loading before proceeding with the placement of a permanent or provisional prosthesis. All components that are used intraorally should be secured to prevent aspiration or swallowing. Distribution of stress is an important consideration. Care should be taken to avoid excessive loads significantly transverse to the implant axis.

## PROSTHETIC COMPONENTS

Preat Prosthetic Components, consisting of abutments, screws, analogs, impression copings, and related restorative accessories, are manufactured from titanium alloy, gold alloy, or polymers. Prosthetic Components are shipped non-sterile (except for Multi-Unit abutments). For product-specific descriptions and sterility information, please contact Preat Corporation.

## Indications for Use

Esthetic Abutments are pre-manufactured components designed as connectors between endosseous dental implants and single unit restorations.

Multi-Unit Abutments are intended to be connectors between endosseous dental implants and multiple implant screw-retained restorations. The 30-degree Multi-Unit abutments must be used within 45 degrees of parallelism. The 17-degree Multi-Unit abutments must be used within 32 degrees of parallelism.

## Impression Procedure

An impression procedure is required whenever a new prosthesis is going to be fabricated. Based on the clinician's preference, open or closed tray impression copings may be used.

**Step 1: Seat the Copings** — Thread and fully seat the closed tray impression coping into the Implant.

If necessary, make adjustments to the tray to allow full seating of the open tray impression copings. Seat the copings and thread the long pins into place. Verify the impression coping is fully seated via x-ray.

**Step 2: Standard impression techniques** are used with the impression copings, recording each implant's position easily and accurately.

**Step 3: Remove the Impression** — Once the impression has fully set, carefully remove the tray from the patient's mouth and verify that all closed tray impression copings have been captured accurately in the impression. Unthread the long guide pins of the open tray impression copings prior to removing the impression tray.

**Step 4: Insert the Analog** — This step can be performed in the clinic or at the dental laboratory. Make sure the closed tray impression coping is free of impression material, and thread an implant analog on to the coping. Prepare the impression to be used to fabricate a stone model.

Thread the implant analog into the long guide pin of the open tray impression coping, assuring full seating. Prepare the impression to be used to fabricate a stone model.

**Step 5: Fabricate the Model** — Use standard laboratory procedure to fabricate a soft tissue model.

## Warnings

Abutments are intended to be used on an individual patient only. The reuse of such device on another patient is not recommended due to the risks of cross-contamination or infection.

## Contraindications

### Esthetic Abutments

- Wall thickness less than 0.5 mm
- Gingival margin diameter less than 0.5 mm wider than implant
- Angle corrections of more than 20 degrees
- Less than 0.5 mm margin height
- Less than 4 mm abutment height

### Multi-Unit Abutments

- Greater than 45 degrees divergence when using 30-degree Multi-Unit abutments
- Greater than 32 degrees divergence when using 17-degree Multi-Unit abutments

## Side Effects

No side effects.

## Sterility

Multi-Unit Abutments are shipped sterile. They should not be resterilized. They are for single use only, prior to the expiration date.

Non-sterile abutments and screws must be cleaned, disinfected, and sterilized prior to clinical use, according to a validated method. The recommended disinfection process is based on ANSI/AAMI 5T79 guidelines, as follows:

- Disinfection: Immerse abutments in disinfectant, rinse with distilled water and dry.
- The recommended sterilization process is based on the ANSI/AAMI/ISO 17665-1 and ANSI/AAMI 5T79 guidelines, as follows:
  - Sterilization: Gravity-fed sterilizers: Autoclave in sterilization pouch for 30 minutes at 121°C (250°F).

NOTE: The validated procedures require the use of FDA-cleared sterilization trays, wraps, biological indicators, chemical indicators, and other sterilization accessories labeled for the sterilization cycle recommended. The healthcare facility should monitor the sterilizer for the facility according to an FDA-recognized sterility assurance standard such as ANSI/AAMI ST79.

## Adverse Effects

The following adverse effects have been observed when using prosthetic components and accessories:

- Components used in the patient's mouth have been aspirated or swallowed.
- The abutment screw has fractured due to application of excessive torque.
- The abutment is not adequately secured due to inadequate application of torque.

## Precautions

Abutments may only be used for their intended purpose in accordance with general rules for prosthetic treatment, occupational safety, and accident prevention. Abutments must only be used for dental procedures with the implant systems they were designed for. If the indications and intended use are not clearly specified, treatment should be suspended until these considerations have been clarified. All components that are used intraorally must be secured to prevent aspiration or swallowing. Prior to placement, ensure that the required components, instruments, and ancillary materials are complete, functional, and available in the correct quantities.

## MRI

Refer to the disinfectant manufacturer's instructions. ANSI/AAMI ST79

Abutments have not been evaluated for safety and compatibility in the MR environment, and have not been tested for heating or migration in the MR environment. They can distort images obtained via magnetic resonance imaging (MRI).

## Dental Implant Compatibility and Recommended Torque Values

Prosthetic Components are generally compatible with the implant systems listed in the table below. The availability of a particular type of prosthetic component varies by implant system. The platform-specific compatibility of each component is indicated on the individual product label.

Implant System and Platform Sizes	Manufacturer's Recommended Torque (Ncm)		
	Titanium Abutment Screw	Multi-Unit Abutment Screw	Multi-Unit Prosthetic Screw
<b>Biomet 3I Certain</b> 3.4 mm, 4.1 mm, 5.0 mm, 6.0 mm	20Ncm	20Ncm	15Ncm
<b>Camlog®Screw-Line</b> 3.3 mm, 3.8 mm, 4.3 mm, 5.0 mm, 6.0 mm	20Ncm	20Ncm	15Ncm
<b>Dentsply Implants Astra Tech Implant System</b> X-Small (2.5 mm) 15 Ncm Small (2.9 mm) 20 Ncm Large (3.9 mm) 25 Ncm	15Ncm 20Ncm 25Ncm	15Ncm 20Ncm 25Ncm	15Ncm 15Ncm 15Ncm
<b>Inclusive Tapered Implant System</b> 3.5 mm, 4.5 mm	35Ncm	30Ncm	15Ncm
<b>Nobel Biocare Brånemark System</b> RP (3.75 mm)	35Ncm	35Ncm	15Ncm
<b>Nobel Biocare NobelActive</b> NP (3.5 mm), RP (4.3 mm)	35Ncm	35Ncm	15Ncm
<b>Nobel Biocare NobelReplace</b> NP (3.5 mm), RP (4.3 mm), WP (5.0 mm), 6.0 (6.0 mm)	35Ncm	35Ncm	15Ncm
<b>Straumann Bone Level</b> NC (3.3 mm), RC (4.1 mm)	35Ncm	35Ncm	15Ncm
<b>Zimmer Dental Screw-Vent®</b> 3.5 mm, 4.5 mm, 5.7 mm	30Ncm	30Ncm	15Ncm

*NOTE: Any screw-retained prosthetic component not listed in the table above should be hand-tightened only.*

ASTRA Tech IMPLANT SYSTEMS is a registered trademark of Dentsply AB LLC. BIOMET 3i™ is a trademark of BIOMET 3i, LLC. Brånemark System™ is a registered trademark of the Nobel Biocare group. CAMLOG™ is a registered trademark of Camlog Biotechnologies AG. CERTAIN™ is a registered trademark of BIOMET 3i, LLC. NobelActive™ is a registered trademark of the Nobel Biocare group. NobelReplace™ is a registered trademark of the Nobel Biocare group. SCREW-VENT™ is a registered trademark of Zimmer Dental Inc. STRAUMANN™ is a registered trademark of Straumann Holding AG.

## INSTRUCTIONS FOR USE — ESTHETIC TITANIUM ABUTMENTS

Esthetic Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to an endosseous implant for retention of a cemented dental prosthesis. They may be indicated for single- and multiple-tooth restorations. Titanium abutments are machined from titanium alloy and attached to the implant fixture with a titanium screw compatible with the restorative instrumentation of the specified implant system.

### Capture Implant Placement

Take an implant-level impression utilizing the preferred technique (open tray, closed tray, or intraoral scan). Submit the impression to the laboratory.

### Laboratory — Fabricate the Restoration

- Follow pouring procedures for the appropriate die stone to produce a working model with implant analog and articulate with a bite registration.
- Select the appropriate Abutment based on the system, platform size, location, and occlusal clearance.
- Seat the abutment completely into the implant analog on the working model, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the emergence profile are esthetically oriented.
- Insert the Titanium Screw (provided) into the abutment's screw access hole and hand-tighten using the appropriate driver.
- Fabricate the restoration using conventional casting or CAD/CAM techniques. Veneer as necessary. If a screw-retained hybrid restoration is indicated, lute the ceramic crown to the titanium abutment.

### Deliver the Final Restoration

- Seat the titanium abutment or screw-retained hybrid restoration completely into the implant, making sure that the anti-rotational features of the connection interface are fully engaged and the contours of the sculpted emergence profile are esthetically oriented.
- Insert the Titanium Screw (provided) into the screw access hole and hand-tighten using the appropriate driver. It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment or hybrid restoration before proceeding.
- Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the abutment or hybrid restoration to the implant manufacturer's recommended torque value.
- Fill the screw access hole with cotton, Teflon tape, gutta-percha, or other suitable material.
- If the restoration is of a screw-retained hybrid design, cover the screw access hole with flowable composite, and cure. Otherwise, follow applicable cementation procedures to affix the definitive restoration to the abutment.

## INSTRUCTIONS FOR USE — MULTI-UNIT ABUTMENTS

Multi-Unit Abutments are prefabricated, screw-retained intraoral abutments intended to be connected directly to endosseous implants in partially or fully edentulous patients for the retention of cast or milled bar overdentures. For implant-supported prostheses, six or more implants are recommended in the maxilla, four or more in the mandible. If clinical conditions dictate fewer implants, an implant-retained, tissue-supported prosthesis is indicated. Multi-unit abutments are machined from titanium alloy, and are available with a variety of collar heights to achieve optimal emergence from shallow or deep gingival wells. Each Multi-Unit Abutment is delivered sterile, suspended in an aseptic vial from a plastic carrier color-coded to indicate the restorative platform of the seated implant.

Straight Multi-Unit abutments lack any anti-rotational features at the implant-abutment interface. The apical portion of a straight Multi-Unit abutment is threaded for integration with the internal cavity of a seated implant. For abutment delivery, the occlusal surface features a male hex head compatible with the Multi-Unit driver recommended by Preat Corporation. Angled Multi-Unit abutments of 17 degrees or 30 degrees enable clinicians and technicians to compensate for the divergence of implants or to otherwise accommodate an angled path of insertion. Angled Multi-Unit abutments feature an anti-rotational connection interface specific to the matching implant platform, and are attached to the implant fixture with an angled Multi-Unit abutment screw compatible with the restorative instrumentation of the specified implant system. Both straight and angled Multi-Unit abutments feature a female connection port at the coronal apex, to allow for the attachment of a screw-retained or fixed-removable dental prosthesis with a Multi-Unit restorative screw.

The axial tilt of an Angled Multi-Unit Abutment (angular divergence from path of insertion) is designed and manufactured to lie along a plane of the implant connection geometry, as opposed to a corner or junction. To maximize the angle-correcting attributes of the Multi-Unit abutment, be sure to rotate the implant upon final seating so that one side of the internal connection geometry (flat or lobe) is oriented to serve as the base of angulation, in accordance with the restorative treatment plan.

### **Place the Multi-Unit Abutment**

- Select the appropriate Multi-Unit Abutment based on platform size, endosseous implant angle and depth of the soft-tissue well.
- Remove the lid from the aseptic vial and retrieve the abutment by lifting the plastic abutment carrier straight out. To maintain the sterility of the Multi-Unit abutment, be careful to handle only by the plastic carrier.
- For Straight Abutments:
  1. Using the plastic carrier, seat the abutment into the implant and hand-tighten. Remove the plastic carrier by pulling the apex of the carrier toward the facial.
- For Angled Abutments:
  1. Using the plastic carrier, seat the abutment into the implant until the anti-rotational features of the connection interface are engaged. Lift and rotate as necessary to orient the angle in the required direction. Hand-tighten the Angled Multi-Unit Abutment Screw using the appropriate driver. Twist the plastic carrier counterclockwise to remove.
  2. NOTE: It is strongly recommended that a radiograph of the connection site be taken to confirm complete seating of the abutment before proceeding.
- Using the appropriate driver in conjunction with a properly metered torque wrench, tighten the Multi-Unit abutment or angled Multi-Unit abutment screw to the implant manufacturer's recommended torque value.

## Passive Temporization of Multi-Unit Abutments

- If the initial stability of the seated implant is insufficient for loading, cover each Multi-Unit Abutment with a Multi-Unit Temporary Healing Cap and hand-tighten with the Prosthetic Screw provided, using the appropriate driver. Do not overtighten.
- Using the patient's existing denture or other prosthesis, relieve the area directly above the placement of each temporary healing cap until the denture rests on the ridge.
- Follow procedures to reline the denture over the temporary healing caps, using soft reline material only. The temporized denture can be used during a healing phase until the implants obtain sufficient load-bearing stability.

## Capture Multi-Unit Abutment Placement

When stability permits, take an abutment-level impression utilizing the preferred technique (open tray, closed tray, or intraoral scan). Submit the impression to the laboratory for the fabrication of a working cast and verification index.

## Denture Protocol

Follow appropriate denture protocol in accordance with the patient-specific treatment plan. When trying in the various setups (e.g., verification index, occlusal rim, wax setup, retention bar), hand-tighten to the Multi-Unit abutments with prosthetic screws, using the appropriate driver. Start from the distal and move forward, alternating between sides of the ridge. Always confirm complete, passive seating, modifying the setup as needed.

## Deliver the Final Restoration

- Remove any temporary prosthesis.
- Confirm that each Multi-Unit abutment is tightened to the implant manufacturer's recommended torque value.

For attachment-retained removable prosthesis:

- Follow procedures to seat the attachment component onto each Multi-Unit abutment. Tighten to the manufacturer's recommended torque value.
- Line the prosthesis onto the attachment components and snap into place. Check comfort and occlusion, and make any necessary adjustments.

For screw-retained, fixed removable prosthesis:

- Line the prosthesis onto the abutments. Beginning with the midmost screw access channel, hand-tighten a Prosthetic Screw into the Multi-Unit abutment. Repeat for each abutment, working outward and alternating left to right.
- Confirm appropriate seating. With the same middle-out, left-to-right technique, tighten each prosthetic screw to 15 Ncm.
- Check comfort and occlusion, and make any necessary adjustments.



- Fill each screw access channel with gutta-percha, silicone, or

For more information, contact Preat at 1-800-232-7732 or visit [preat.com](http://preat.com)